

AW Technologies Press Release / 09.02.2026


# AW Technologies and Diamed Medizintechnik announces they have entered into an exclusive distribution partnership for TrachFlush™ in Germany

Nørresundby, Denmark / Cologne, Germany – AW Technologies, a company focused on developing and bringing to market disruptive and innovative airway solutions for mechanically ventilated patients, and Diamed, a company focused on bringing innovative extracorporeal blood purification and airway management to the German market, announced today they have entered into an exclusive distribution agreement for the TrachFlush™ device in Germany.

Invented by Alberto Zanella, an Italian intensivist, and developed by AW Technologies, TrachFlush™ is a medical device that maintains the tracheal tube cuff's pressure, preventing subglottic secretion microaspiration, and removes the secretion below the cuff and above (subglottic) the cuff non-invasively, lightening the workload for healthcare professionals and reduced patient discomfort.



**The TrachFlush™ Device**



**Illustration of removing secretions with TrachFlush**

When the Flush button is pushed, TrachFlush automatically synchronizes with the ventilator, and deflates the cuff during the inspiratory cycle, utilizing the airflow of pressure from the lungs and the ventilator to “push” the secretions from above and below the cuff up into the mouth of the patient. TrachFlush re-inflates the cuff again prior to the expiratory cycle, avoiding PEEP loss.

**“With AW Technologies, we are achieving another milestone in our mission to advance innovative therapy concepts for intensive care medicine. Together with our agile partner, we aim to help improve patient care quality and ease the clinical environment through new, practical solutions”**

MR. JOACHIM SEYFANG

Chief Executive Officer of Diamed Medizintechnik GmbH

**“We are excited to partner with DIAMED to bring TrachFlush to the German market. Together, we aim to expand access to effective, safe, and high quality care through shared expertise and commitment to patient well-being.**

MR. ADAM HANSEN

Chief Executive Officer of AW Technologies ApS

### **About Diamed Medizintechnik GmbH**

Diamed Medizintechnik GmbH has been established in intensive care medicine for 50 years and provides solutions for extracorporeal blood purification, as well as newly for secretion removal in invasive ventilation.

For more information about DIAMED, please visit <https://diamed.de/home.html>.

### **About AW Technologies ApS**

AW Technologies – a Danish medical device company, established in late November 2018, focused on developing and bringing to market disruptive and innovative airway solutions for mechanically ventilated patients, including the primary product; TrachFlush. For more information, please visit [www.trachflush.com](http://www.trachflush.com) or [www.awtechnologies.dk](http://www.awtechnologies.dk)

### **Contact**

AW Technologies ApS  
Amalienborgvej 57  
9400, Nørresundby, Denmark  
+45 7777 4048  
info@awtechnologies.dk  
[www.awtechnologies.dk](http://www.awtechnologies.dk)  
[www.trachflush.com](http://www.trachflush.com)